

# BeneFusion eSP/eVP/eDS Vet

Veterinary Syringe Pump / Infusion Pump / Infusion Supervision System

## Efficiency in every droplet



LinkedIn Facebook Instagram

Search Mindray Animal Medical

[www.mindrayanimal.com/na/](http://www.mindrayanimal.com/na/)


P/N:ENG-BeneFusion eSPeVPeDS Vet-210285X6P-20230213

© 2023 Shenzhen Mindray Animal Medical Technology Co., LTD. All rights reserved.

## Efficiency in Workflow


### SmartRapid™

SmartRapid™ ensures timely infusion by significantly shortening the start-up time, from turning on the pump to delivering the first drop of medication.




## Intuitive interaction

- 3.5" Colored capacitive touchscreen brings users an excellent experience with smooth operation.
- Intuitive user interface enables quick programming of key parameters, making workflow process much easier.



## Flexible infusion adjustment


Easily adjusting the infusion therapy within 3 steps.



## Efficiency in Safety

### Precise infusion performance

- High accuracy: eSP ±1.8% / eVP ±5%
- Long-hour accuracy assured



### SmartAIR™

With SmartAIR™, BeneFusion eVP Vet significantly enhances the safety level of IV administration.

- Dual ultrasound sensor to detect the air bubbles more precisely, avoiding missing or false air-in-line alarm.
- 15µl air bubble detection size on BeneFusion eVP Vet, ensuring patient safety throughout the infusion, even for neonates.



## Color coding

- The color coding of drug name assists users to easily select and verify the correct drug.



## Quick problem solving

### Dynamic Pressure System

Speedometer style indicator with numerical pressure to monitor the in-line pressure trend just with a glance.




Alert with a visual message about the possible occlusion before interrupting the infusion.

Auto-attempting mechanism enables infusion resume as early as possible to ensure continuous infusion.

### Intuitive Alarm System

Instructional animation pops up to guide users to quickly solve the problem.



## Efficiency in Application

### All in one

BeneFusion eSP and eVP Vet satisfy various infusion purposes by combining all functions together.



## Integrated central monitoring

BeneVision CMS Vet offers one-stop monitoring of all patients' vital signs and infusion treatment details, providing comprehensive information for the veterinarian to improve the quality of patient care.



Easy management with multi-beds



Comprehensive data for single bed