

youlan

YL VET DENTAL UNIT ProDx Series

Enjoy Health, Enjoy life



youlan
BETTER CARE FOR ANIMALS

SHENZHEN YOULAN MEDICAL TECHNOLOGY CO., Ltd.

✉ Email: sales@ylvet.com

📍 Address: Fifth Floor, Building 4, Baoyunda Logistics Center, Fuhua Community,
Xixiang Street, Bao'an District, Shenzhen, Guangdong Province, China

☎ 0755-23054282

🌐 www.ylvet.com



YL VET DENTAL UNIT

ProDx Series

ProDx Series workstation can serve as a comprehensive dental work platform for tooth polishing, restoration, cutting and root canal treatment. It can also be used for animal tooth beautification, tooth removal, oral tumor cutting, etc., bringing more help to animal welfare. It can also improve the veterinarian's handling of animal oral problems more efficiently. With a higher hardware selection standard, open design, and customized hardware configuration, this model brings the ultimate experience to veterinary practice.

Multifunctional dental equipment



BETTER CARE FOR ANIMALS 🦷



Smart & Efficient

Design for Animals

YL VET DENTAL UNIT

ProDx Series



Smart

Area size 44*36.5cm, saving space; Height is adjustable, between 90cm and 105cm, suitable for the diverse practice requirements of veterinarians.



High efficiency and low noise

High-quality oil-free air compressor, lower noise. Water pressure flow can be controlled simply and efficiently

BETTER CARE FOR ANIMALS 

Product Info.

Description	
Compressor	1pcs
Power Line	1pcs
High Speed Handpiece	1pcs
Low Speed Handpiece	1pcs
Scaler Handpiece	1pcs
3-Way Air/Water Syringe	1pcs
Burs	1 Set
Disposable Polishing Head	1 Set
Scaling Tips	1 Set (5pcs)

Technical Parameters

Item	Specification
Compressor	Oil/Oil-free Compressor
Power	2.6-2.9A, 220/230V, 50/60Hz
Output Pressure	0.45~0.5MPa
Air Tank	8/20L(Customized)
Sound Level	<65dB (Oil-Free); <50dB(Oil-Compressor)
High Speed Handpiece	300,000~340,000rpm
Low Speed Handpiece	1800~20000rpm
Ultrasonic Scaler	With White LED light 28kHz±3kHz
Dimension	L44*W36.5*H76-116cm
Weight	About 40kg